



Airways - Oropharyngeal



COLOR-CODED GUEDEL

- Color coded, disposable
- Firm airway
- Built-in bite block
- Latex free - polyethylene plastic
- **Individually wrapped, non-sterile**

COLOR-CODED GUEDEL AIRWAY SET

- Ideal for crash carts
- Individually wrapped, non-sterile

Designed and manufactured in accordance with ANSI standard for oropharyngeal airways.



COLOR-CODED PVC GUEDEL

- Soft, clear PVC plastic
- **Individually wrapped, non-sterile**
- Latex free

COLOR-CODED PVC GUEDEL AIRWAY SET

- Ideal for crash carts
- **Individually wrapped, non-sterile**

CAT#	SIZE	DESCRIPTION	COLOR	BOX
1-1504-40	00	Newborn, 40mm	Pink	10
1-1504-50	0	Infant, 50mm	Blue	10
1-1504-60	1	Small Child, 60mm	Black	10
1-1504-70	2	Child, 70mm	White	10
1-1504-80	3	Small Adult, 80mm	Green	10
1-1504-90	4	Med. Adult, 90mm	Yellow	10
1-1504-99	5	Large Adult, 100mm	Red	10
1-1504-11	6	X-Large Adult, 110mm	Orange	10
1-1505-40	00	Newborn, 40mm	Pink	50
1-1505-50	0	Infant, 50mm	Blue	50
1-1505-60	1	Small Child, 60mm	Black	50
1-1505-70	2	Child, 70mm	White	50
1-1505-80	3	Small Adult, 80mm	Green	50
1-1505-90	4	Med. Adult, 90mm	Yellow	50
1-1505-99	5	Large Adult, 100mm	Red	50
1-1505-11	6	X-Large Adult, 110mm	Orange	50

COLOR CODED GUEDEL AIRWAY SET	BOX
1-1504-00	One of each size: 40, 50, 60, 70, 80, 90, 100, 110mm
	1 Set

CAT#	SIZE	DESCRIPTION	COLOR	BOX
1-5010-30	000	Premature, 30mm	Clear	10
1-5010-20	00	Newborn, 40mm	Pink	10
1-5010-10	0	Infant, 50mm	Blue	10
1-5010-01	1	Small Child, 60mm	Black	10
1-5010-02	2	Child, 70mm	White	10
1-5010-03	3	Small Adult, 80mm	Green	10
1-5010-04	4	Med. Adult, 90mm	Yellow	10
1-5010-05	5	Large Adult, 100mm	Red	10
1-5010-06	6	X-Large Adult, 110mm	Orange	10
1-1521-30	000	Premature, 30mm	Clear	50
1-1521-20	00	Newborn, 40mm	Pink	50
1-1521-10	0	Infant, 50mm	Blue	50
1-1521-01	1	Small Child, 60mm	Black	50
1-1521-02	2	Child, 70mm	White	50
1-1521-03	3	Small Adult, 80mm	Green	50
1-1521-04	4	Med. Adult, 90mm	Yellow	50
1-1521-05	5	Large Adult, 100mm	Red	50
1-1521-06	6	X-Large Adult, 110mm	Orange	50

COLOR CODED PVC GUEDEL AIRWAY SET	BOX
1-5010-98	One of each size: 30, 40, 50, 60, 70, 80, 90, 100 & 110mm
	1 Set




Airways - Oralpharyngeal



TRADITIONAL GUEDEL AIRWAYS

- Latex free - polyethylene plastic
- **Not individually wrapped**



WILLIAMS AIRWAY INTUBATOR

The pink airway intubator is indicated for use as:

- An oropharyngeal airway
- A means of intubating the trachea
- A guide for fiber optic bronchoscopy placement
- Disposable and latex free
- Built-in bite block



COLOR CODED BERMAN AIRWAYS

- Color coded for easy identification
- Latex free - polyethylene plastic
- Vented
- **Disposable and individually wrapped**

CAT#	SIZE	DESCRIPTION	COLOR	BOX
1-1501-50	0	Infant, 50mm	Blue	10
1-1501-60	1	Small Child, 60mm	Black	10
1-1501-70	2	Child, 70mm	White	10
1-1501-80	3	Small Adult, 80mm	Green	10
1-1501-90	4	Medium Adult, 90mm	Yellow	10
1-1501-99	5	Large Adult, 100mm	Red	10
1-1500-50	0	Infant, 50mm	Blue	50
1-1500-60	1	Small Child, 60mm	Black	50
1-1500-70	2	Child, 70mm	White	50
1-1500-80	3	Small Adult, 80mm	Green	50
1-1500-90	4	Medium Adult, 90mm	Yellow	50
1-1500-99	5	Large Adult, 100mm	Red	50

TRADITIONAL GUEDEL DISPOSABLE AIRWAY SET	BOX
1-1501-00 Set of one of each size: 50, 60, 70, 80, 90, & 100mm	1 Set

CAT#	SIZE	DESCRIPTION	BOX
1-1503-09	9cm	Adult Female	10
1-1503-10	10cm	Adult Male	10

CAT#	SIZE	DESCRIPTION	COLOR	BOX
1-1508-40	00	Newborn, 40mm	Pink	10
1-1508-50	0	Infant, 50mm	Blue	10
1-1508-60	1	Child, 60mm	Black	10
1-1508-70	2	Small Child, 70mm	White	10
1-1508-80	3	Small Adult, 80mm	Green	10
1-1508-90	4	Medium Adult, 90mm	Yellow	10
1-1508-99	5	Large Adult, 100mm	Red	10
1-1508-11	6	X-Large Adult, 110mm	Orange	10
1-1506-40	00	Newborn, 40mm	Pink	50
1-1506-50	0	Infant, 50mm	Blue	50
1-1506-60	1	Child, 60mm	Black	50
1-1506-70	2	Small Child, 70mm	White	50
1-1506-80	3	Small Adult, 80mm	Green	50
1-1506-90	4	Medium Adult, 90mm	Yellow	50
1-1506-99	5	Large Adult, 100mm	Red	50
1-1506-11	6	X-Large Adult, 110mm	Orange	50

COLOR CODED BERMAN DISPOSABLE AIRWAY SET	BOX
1-1508-00 Set of one of each size: 40, 50, 60, 70, 80, 90, 100 & 110mm	10 Sets



Airways - Nasopharyngeal



Latex Free



ROBERTAZZI (TRUMPET) STYLE

Latex free nasopharyngeal airways have the same characteristics as the popular "Blue Latex" airway, in softness and flexibility.

- Latex free
- Soft - no damage to nasal passage
- Flexible
- Resists collapse or kinking
- Resilient
- Nine sizes available: 20-36FR
- Sterile & disposable

CAT#	SIZE	DESCRIPTION	BOX
1-5075-20	20FR	Robertazzi Style Airway	10
1-5075-22	22FR	Robertazzi Style Airway	10
1-5075-24	24FR	Robertazzi Style Airway	10
1-5075-26	26FR	Robertazzi Style Airway	10
1-5075-28	28FR	Robertazzi Style Airway	10
1-5075-30	30FR	Robertazzi Style Airway	10
1-5075-32	32FR	Robertazzi Style Airway	10
1-5075-34	34FR	Robertazzi Style Airway	10
1-5075-36	36FR	Robertazzi Style Airway	10

STERILE - DISPOSABLE - LATEX FREE NASOPHARYNGEAL KITS

1-5075-04	Sizes 20, 24, 28, 32FR w/sterile lubricating jelly	4/pack
1-5075-06	Sizes 22, 24, 26, 28, 30, 32FR w/sterile lubricating jelly	6/pack
1-5075-09	Sizes 20, 22, 24, 26, 28, 30, 32, 34, 36 FR w/sterile lubricating jelly	9/pack



ADJUSTABLE FLANGE

- Soft moveable flange
- Thirteen sizes 12-36FR
- Latex free
- Sterile & disposable

CAT#	SIZE	DESCRIPTION	BOX
1-5072-12	12FR	Adjustable Flange Airway	10
1-5072-14	14FR	Adjustable Flange Airway	10
1-5072-16	16FR	Adjustable Flange Airway	10
1-5072-18	18FR	Adjustable Flange Airway	10
1-5072-20	20FR	Adjustable Flange Airway	10
1-5072-22	22FR	Adjustable Flange Airway	10
1-5072-24	24FR	Adjustable Flange Airway	10
1-5072-26	26FR	Adjustable Flange Airway	10
1-5072-28	28FR	Adjustable Flange Airway	10
1-5072-30	30FR	Adjustable Flange Airway	10
1-5072-32	32FR	Adjustable Flange Airway	10
1-5072-34	34FR	Adjustable Flange Airway	10
1-5072-36	36FR	Adjustable Flange Airway	10



Airways - Nasopharyngeal



Original Style "Blue" Robertazzi

"BLUE" LATEX ROBERTAZZI (TRUMPET)

FEATURES

- Sterile and flexible
- Smooth surface for easy insertion
- Soft - to reduce trauma to nasal passage
- Resists collapse or kinking
- Process to reduce latex protein




CAT#	SIZE	DESCRIPTION	BOX
1-5076-20	20FR	Robertazzi "Blue" Latex	10
1-5076-22	22FR	Robertazzi "Blue" Latex	10
1-5076-24	24FR	Robertazzi "Blue" Latex	10
1-5076-26	26FR	Robertazzi "Blue" Latex	10
1-5076-28	28FR	Robertazzi "Blue" Latex	10
1-5076-30	30FR	Robertazzi "Blue" Latex	10
1-5076-32	32FR	Robertazzi "Blue" Latex	10
1-5076-34	34FR	Robertazzi "Blue" Latex	10
1-5076-36	36FR	Robertazzi "Blue" Latex	10

STERILE - "BLUE" LATEX NASOPHARYNGEAL KITS

1-5076-04	Sizes 20, 24, 28, 32FR w/sterile lubricating jelly	4/pack
1-5076-06	Sizes 22, 24, 26, 28, 30, 32FR w/sterile lubricating jelly	6/pack
1-5076-09	Sizes 20, 22, 24, 26, 28, 30, 32, 34, 36 FR w/sterile lubricating jelly	9/pack

PVC AIRWAY

- Sterile and Latex free
- Beveled Tip
- Single use
- 15mm Fitting



CAT#	SIZE	DESCRIPTION	BOX
1-5073-20	5mm/20FR	PVC Naso	10
1-5073-24	6mm/24FR	PVC Naso	10
1-5073-28	7mm/28FR	PVC Naso	10
1-5073-32	8mm/32FR	PVC Naso	10
1-5073-36	9mm/36FR	PVC Naso	10

PVC NASOPHARYNGEAL KITS


1-5073-05	Sizes 20, 24, 28, 32, 36FR w/sterile lubricating jelly	5/pack
-----------	---	--------



Airways - Oral and Nasal



Endotracheal Tubes



MURPHY UN-CUFFED				
CAT#	I.D.	FRENCH	LENGTH	BOX
1-7330-20	2.0mm	8FR	140mm	10
1-7330-25	2.5mm	10FR	145mm	10
1-7330-30	3.0mm	12FR	165mm	10
1-7330-35	3.5mm	14FR	185mm	10
1-7330-40	4.0mm	16FR	210mm	10
1-7330-45	4.5mm	18FR	225mm	10
1-7330-50	5.0mm	20FR	245mm	10
1-7330-55	5.5mm	22FR	275mm	10

MURPHY CUFFED				
CAT#	I.D.	FRENCH	LENGTH	BOX
1-7333-30	3.0mm	12FR	165mm	10
1-7333-35	3.5mm	14FR	185mm	10
1-7333-40	4.0mm	16FR	210mm	10
1-7333-45	4.5mm	18FR	225mm	10
1-7333-50	5.0mm	20FR	245mm	10
1-7333-55	5.5mm	22FR	275mm	10
1-7333-60	6.0mm	24FR	285mm	10
1-7333-65	6.5mm	26FR	295mm	10
1-7333-70	7.0mm	28FR	310mm	10
1-7333-75	7.5mm	30FR	315mm	10
1-7333-80	8.0mm	32FR	330mm	10
1-7333-85	8.5mm	34FR	330mm	10
1-7333-90	9.0mm	36FR	330mm	10
1-7333-95	9.5mm	38FR	330mm	10
1-7333-10	10.0mm	40FR	330mm	10

The cuffed tubes are manufactured with a large barrel type cuff to increase surface contact for a gentle seal during tracheal expansion and contraction. All endotracheal tubes have a smooth Murphy Eye to provide an alternate ventilation port if tube tip becomes occluded.

- For Oral and Nasal Use
- Available Cuffed or Un-Cuffed
- Most Economical - Sterile & Disposable
- Radio Opaque Strip Embedded
- Smooth Tube Tip to Minimize Trauma During Intubation
- High Volume, Low Pressure Cuff
- Standard 15mm ISO Connector Supplied with Each Tube
- Latex Free - Tubes are Constructed of Clear, Flexible PVC
- Bold, Easy-to-read Depth Marks for Placement in Low Light Conditions
- Valve with Sensitive Pilot Balloon Indicating Status of Cuff




Airways - Oral and Nasal



Endotracheal Tubes With Stylets

CAT#	I.D.	FRENCH	LENGTH	BOX
1-7343-50	5.0mm	20FR	245mm	10
1-7343-55	5.5mm	22FR	275mm	10
1-7343-60	6.0mm	24FR	285mm	10
1-7343-65	6.5mm	26FR	295mm	10
1-7343-70	7.0mm	28FR	310mm	10
1-7343-75	7.5mm	30FR	315mm	10
1-7343-80	8.0mm	32FR	330mm	10
1-7343-85	8.5mm	34FR	330mm	10
1-7343-90	9.0mm	36FR	330mm	10



ENDOTRACHEAL TUBES WITH STYLETS

- Convenient - Tube Supplied with Stylet Inside
- Surface of Stylet Treated to Prevent Friction While Sliding In and Out of Tube
- Most Economical - Sterile, Disposable & Latex Free
- Radio Opaque Strip Embedded
- Smooth Tube Tip to Minimize Trauma During Intubation
- High Volume, Low Pressure Cuff

CUFF INFLATION BULB

- Quickly Inflates and Deflates Endotracheal Tube Cuffs
- Safe - No Syringes
- Latex Free
- Fits Easily In Pocket - One Hand Operation



CAT#	DESCRIPTION	BOX
4-0101-00	Cuff Inflation Bulb (bulb only)	10



Also Inflates and Deflates Mask Cushions