




# Airways - Oral and Nasal



## Endotracheal Tubes



MURPHY UN-CUFFED				
CAT#	I.D.	FRENCH	LENGTH	BOX
1-7330-20	2.0mm	8FR	140mm	10
1-7330-25	2.5mm	10FR	145mm	10
1-7330-30	3.0mm	12FR	165mm	10
1-7330-35	3.5mm	14FR	185mm	10
1-7330-40	4.0mm	16FR	210mm	10
1-7330-45	4.5mm	18FR	225mm	10
1-7330-50	5.0mm	20FR	245mm	10
1-7330-55	5.5mm	22FR	275mm	10

MURPHY CUFFED				
CAT#	I.D.	FRENCH	LENGTH	BOX
1-7333-30	3.0mm	12FR	165mm	10
1-7333-35	3.5mm	14FR	185mm	10
1-7333-40	4.0mm	16FR	210mm	10
1-7333-45	4.5mm	18FR	225mm	10
1-7333-50	5.0mm	20FR	245mm	10
1-7333-55	5.5mm	22FR	275mm	10
1-7333-60	6.0mm	24FR	285mm	10
1-7333-65	6.5mm	26FR	295mm	10
1-7333-70	7.0mm	28FR	310mm	10
1-7333-75	7.5mm	30FR	315mm	10
1-7333-80	8.0mm	32FR	330mm	10
1-7333-85	8.5mm	34FR	330mm	10
1-7333-90	9.0mm	36FR	330mm	10
1-7333-95	9.5mm	38FR	330mm	10
1-7333-10	10.0mm	40FR	330mm	10

The cuffed tubes are manufactured with a large barrel type cuff to increase surface contact for a gentle seal during tracheal expansion and contraction. All endotracheal tubes have a smooth Murphy Eye to provide an alternate ventilation port if tube tip becomes occluded.

- For Oral and Nasal Use
- Available Cuffed or Un-Cuffed
- Most Economical - Sterile & Disposable
- Radio Opaque Strip Embedded
- Smooth Tube Tip to Minimize Trauma During Intubation
- High Volume, Low Pressure Cuff
- Standard 15mm ISO Connector Supplied with Each Tube
- Latex Free - Tubes are Constructed of Clear, Flexible PVC
- Bold, Easy-to-read Depth Marks for Placement in Low Light Conditions
- Valve with Sensitive Pilot Balloon Indicating Status of Cuff




# Airways - Oral and Nasal



## Endotracheal Tubes With Stylets

CAT#	I.D.	FRENCH	LENGTH	BOX
1-7343-50	5.0mm	20FR	245mm	10
1-7343-55	5.5mm	22FR	275mm	10
1-7343-60	6.0mm	24FR	285mm	10
1-7343-65	6.5mm	26FR	295mm	10
1-7343-70	7.0mm	28FR	310mm	10
1-7343-75	7.5mm	30FR	315mm	10
1-7343-80	8.0mm	32FR	330mm	10
1-7343-85	8.5mm	34FR	330mm	10
1-7343-90	9.0mm	36FR	330mm	10



### ENDOTRACHEAL TUBES WITH STYLETS

- Convenient - Tube Supplied with Stylet Inside
- Surface of Stylet Treated to Prevent Friction While Sliding In and Out of Tube
- Most Economical - Sterile, Disposable & Latex Free
- Radio Opaque Strip Embedded
- Smooth Tube Tip to Minimize Trauma During Intubation
- High Volume, Low Pressure Cuff

### CUFF INFLATION BULB

- Quickly Inflates and Deflates Endotracheal Tube Cuffs
- Safe - No Syringes
- Latex Free
- Fits Easily In Pocket - One Hand Operation



CAT#	DESCRIPTION	BOX
4-0101-00	Cuff Inflation Bulb (bulb only)	10



Also Inflates and Deflates Mask Cushions



# Airways



## Tracheostomy Tubes

- Obturator and Tube Clearly Marked
- 15mm Connector Swivel Adapter
- Securely Attached Neck Band Retainer
- X-ray Opaque
- Low Pressure Cuff
- Comfortable Cushion Neck Band
- Made of Soft PVC
- Latex Free - Single Use
- Sterile



CAT#	DESCRIPTION	SIZE	BOX
1-7390-30	Trach Tube Uncuffed	3.0mm	10
1-7390-35	Trach Tube Uncuffed	3.5mm	10
1-7390-40	Trach Tube Uncuffed	4.0mm	10
1-7390-45	Trach Tube Uncuffed	4.5mm	10
1-7390-50	Trach Tube Uncuffed	5.0mm	10



CAT#	DESCRIPTION	SIZE	BOX
1-7393-50	Trach Tube Cuffed	5.0mm	10
1-7393-55	Trach Tube Cuffed	5.5mm	10
1-7393-60	Trach Tube Cuffed	6.0mm	10
1-7393-65	Trach Tube Cuffed	6.5mm	10
1-7393-70	Trach Tube Cuffed	7.0mm	10
1-7393-75	Trach Tube Cuffed	7.5mm	10
1-7393-80	Trach Tube Cuffed	8.0mm	10
1-7393-85	Trach Tube Cuffed	8.5mm	10
1-7393-90	Trach Tube Cuffed	9.0mm	10
1-7393-95	Trach Tube Cuffed	9.5mm	10
1-7393-10	Trach Tube Cuffed	10.0mm	10