



# Urology Products



## Latex Silicone Coated Foley Catheter


- Smooth tapered tip facilitates easy insertion into urethra
- Drainage eyes are accurately formed to permit effective drainage
- Symmetrical balloon expands equally in all directions and efficiently retains bladder
- Color-coded sleeves for easy and rapid size identification
- Plastic valve permits simple inflation and deflation - more efficient than rubber valves

CAT#	DESCRIPTION	SIZE	BOX
7-6503-08	Foley Pedi Latex 2 Way	3cc, 8FR	10
7-6503-10	Foley Pedi Latex 2 Way	3cc, 10FR	10

CAT#	DESCRIPTION	SIZE	BOX
7-6505-12	Foley Latex 2 Way	5cc, 12FR	10
7-6505-14	Foley Latex 2 Way	5cc, 14FR	10
7-6505-16	Foley Latex 2 Way	5cc, 16FR	10
7-6505-18	Foley Latex 2 Way	5cc, 18FR	10
7-6505-20	Foley Latex 2 Way	5cc, 20FR	10
7-6505-22	Foley Latex 2 Way	5cc, 22FR	10
7-6505-24	Foley Latex 2 Way	5cc, 24FR	10
7-6505-26	Foley Latex 2 Way	5cc, 26FR	10
7-6505-28	Foley Latex 2 Way	5cc, 28FR	10
7-6505-30	Foley Latex 2 Way	5cc, 30FR	10

7-6507-12	Foley Latex 2 Way	30cc, 12FR	10
7-6507-14	Foley Latex 2 Way	30cc, 14FR	10
7-6507-16	Foley Latex 2 Way	30cc, 16FR	10
7-6507-18	Foley Latex 2 Way	30cc, 18FR	10
7-6507-20	Foley Latex 2 Way	30cc, 20FR	10
7-6507-22	Foley Latex 2 Way	30cc, 22FR	10
7-6507-24	Foley Latex 2 Way	30cc, 24FR	10
7-6507-26	Foley Latex 2 Way	30cc, 26FR	10
7-6507-28	Foley Latex 2 Way	30cc, 28FR	10
7-6507-30	Foley Latex 2 Way	30cc, 30FR	10

CAT#	DESCRIPTION	SIZE	BOX
7-6509-18	Foley Latex 3 Way	30cc, 18FR	10
7-6509-20	Foley Latex 3 Way	30cc, 20FR	10
7-6509-22	Foley Latex 3 Way	30cc, 22FR	10
7-6509-24	Foley Latex 3 Way	30cc, 24FR	10
7-6509-26	Foley Latex 3 Way	30cc, 26FR	10
7-6509-28	Foley Latex 3 Way	30cc, 28FR	10





# Urology Products



## 100% Silicone Foley Catheter


- Silicone - Inert material preferred for medical devices, less urethral irritation and stricture versus latex material
- Less encrustations compared to latex catheters. 100% latex-free.
- Smooth tapered tip facilitates easy insertion into urethra
- Drainage eyes are accurately formed to permit effective drainage
- Symmetrical balloon expands equally in all directions and efficiently retains bladder
- Color-coded sleeves for easy and rapid size identification
- Plastic Bespak spring loaded valve permits simple inflation and deflation - more efficient than rubber valves

CAT#	DESCRIPTION	SIZE	BOX
7-6513-08	Pedi Silicone 2 Way	3cc, 8FR	10
7-6513-10	Pedi Silicone 2 Way	3cc, 10FR	10

CAT#	DESCRIPTION	SIZE	BOX
7-6515-12	Foley Silicone 2 Way	5cc, 12FR	10
7-6515-14	Foley Silicone 2 Way	5cc, 14FR	10
7-6515-16	Foley Silicone 2 Way	5cc, 16FR	10
7-6515-18	Foley Silicone 2 Way	5cc, 18FR	10
7-6515-20	Foley Silicone 2 Way	5cc, 20FR	10
7-6515-22	Foley Silicone 2 Way	5cc, 22FR	10
7-6515-24	Foley Silicone 2 Way	5cc, 24FR	10
7-6515-26	Foley Silicone 2 Way	5cc, 26FR	10
7-6517-16	Foley Silicone 2 Way	30cc, 16FR	10
7-6517-18	Foley Silicone 2 Way	30cc, 18FR	10
7-6517-20	Foley Silicone 2 Way	30cc, 20FR	10
7-6517-22	Foley Silicone 2 Way	30cc, 22FR	10
7-6517-24	Foley Silicone 2 Way	30cc, 24FR	10
7-6517-26	Foley Silicone 2 Way	30cc, 26FR	10

CAT#	DESCRIPTION	SIZE	BOX
7-6519-16	Foley Silicone 3 Way	30cc, 16FR	10
7-6519-18	Foley Silicone 3 Way	30cc, 18FR	10
7-6519-20	Foley Silicone 3 Way	30cc, 20FR	10
7-6519-22	Foley Silicone 3 Way	30cc, 22FR	10
7-6519-24	Foley Silicone 3 Way	30cc, 24FR	10
7-6519-26	Foley Silicone 3 Way	30cc, 26FR	10

CAT#	DESCRIPTION	SIZE	BOX
7-66614-02	Irrigating Bulb Syringe With Tip	60cc	50





# Urology Products



## Urinary Drainage Bags

- Single hook bed hanger
- Anti-reflux feature
- Vented bag
- Sampling port & needleless access
- 40" clear inlet PVC tubing
- T-tap outlet device
- Sterile & latex free



CAT#	DESCRIPTION	SIZE	BOX
7-6610-00	Bedside Economy Bag	2000ml (64oz.)	20

- Dual hook bed hanger
- String tie for added security
- Drip chamber w/anti-reflux valve
- Thicker more durable vented bag
- Sampling port with needleless access
- 48" clear inlet PVC tubing & sheet clamp
- T-tap outlet device
- Sterile & latex free



CAT#	DESCRIPTION	SIZE	BOX
7-6612-00	Bedside Premium Bag	2000ml (64oz.)	10

- Velcro bed hanger
- Drip chamber with anti-reflux valve
- Vented bag
- 48" clear inlet PVC tubing
- Sampling port with needleless access
- T-tap outlet device
- Sterile & latex free



CAT#	DESCRIPTION	SIZE	BOX
7-6610-50	Metered Bedside Bag 3-Chamber 400ml meter	2000ml (64oz.)	20

- Anti-reflux valve & twist outlet
- Two elastic leg straps are pre installed on the bag
- Two buttons provided to secure straps
- Sterile & Latex free



CAT#	DESCRIPTION	SIZE	BOX
7-6600-19	Leg Bag Medium	500ml (19oz.)	48
7-6600-32	Leg Bag Large	1000ml (32oz.)	48